


euca PW - PTA Balloon Catheter-OTW



flexible tip
reduced
crossing profile
low balloon profile

new balloon folding

Dual lumen shaft technology

Laser rounded tip

Controlled compliance

Excellent refolding

Rapid deflation

New Generation
PTA by eucatech
0.035" Guide wire

0.035" Guide wire

Balloon diameter from 4 to 12 mm

Balloon length from 20 to 200 mm

Sheath compatibility 5 to 8F

36 months shelf life

up to 200 mm
balloon length

Long balloons for optimal lesion treatment.



Quality Made in Germany

euca PW - PTA Balloon Catheter-OTW

Catheter Details	
Usable Length	85 cm / 135 cm
Balloon Material	Proprietary polyamide
Shaft Material	Polyamide, double-lumen tube
X-Ray Markers	platinum/iridium
Max. recommended Guide Wire	0.035" (0.89 mm)

Order Information					
Balloon Ø	Balloon length	RBP	Rec. ID Introducer Sheath	Usable Catheter Length	
[mm]	[mm]	[atm]	[F/mm]	85 cm	135 cm
4	20	14	5/1.65	PW 04/020-085	PW 04/020-135
	30			PW 04/030-085	PW 04/030-135
	40			PW 04/040-085	PW 04/040-135
	60			PW 04/060-085	PW 04/060-135
	80			PW 04/080-085	PW 04/080-135
	100			PW 04/100-085	PW 04/100-135
	120			PW 04/120-085	PW 04/120-135
	150			PW 04/150-085	PW 04/150-135
	170			PW 04/170-085	PW 04/170-135
	200			PW 04/200-085	PW 04/200-135
5	20	13	5/1.65	PW 05/020-085	PW 05/020-135
	30			PW 05/030-085	PW 05/030-135
	40			PW 05/040-085	PW 05/040-135
	60			PW 05/060-085	PW 05/060-135
	80			PW 05/080-085	PW 05/080-135
	100	6/1.98	PW 05/100-085	PW 05/100-135	
	120		PW 05/120-085	PW 05/120-135	
	150		PW 05/150-085	PW 05/150-135	
	170		PW 05/170-085	PW 05/170-135	
	200		PW 05/200-085	PW 05/200-135	
6	20	12	6/1.98	PW 06/020-085	PW 06/020-135
	30			PW 06/030-085	PW 06/030-135
	40			PW 06/040-085	PW 06/040-135
	60			PW 06/060-085	PW 06/060-135
	80			PW 06/080-085	PW 06/080-135
	100			PW 06/100-085	PW 06/100-135
	120			PW 06/120-085	PW 06/120-135
	150			PW 06/150-085	PW 06/150-135
	170			PW 06/170-085	PW 06/170-135
	200			PW 06/200-085	PW 06/200-135
7	20	11	6/1.98	PW 07/020-085	PW 07/020-135
	30			PW 07/030-085	PW 07/030-135
	40			PW 07/040-085	PW 07/040-135
	60			PW 07/060-085	PW 07/060-135
	80			PW 07/080-085	PW 07/080-135
	100	7/2.33	PW 07/100-085	PW 07/100-135	
	120		PW 07/120-085	PW 07/120-135	
	150		PW 07/150-085	PW 07/150-135	
	170		PW 07/170-085	PW 07/170-135	
	200		PW 07/200-085	PW 07/200-135	
8	20	11	7/2.33	PW 08/020-085	PW 08/020-135
	30			PW 08/030-085	PW 08/030-135
	40			PW 08/040-085	PW 08/040-135
	60			PW 08/060-085	PW 08/060-135
	80			PW 08/080-085	PW 08/080-135
9	20	11	7/2.33	PW 09/020-085	PW 09/020-135
	30			PW 09/030-085	PW 09/030-135
	40			PW 09/040-085	PW 09/040-135
	60			PW 09/060-085	PW 09/060-135
10	20	10	7/2.33	PW 10/020-085	PW 10/020-135
	30			PW 10/030-085	PW 10/030-135
	40			PW 10/040-085	PW 10/040-135
	60			PW 10/060-085	PW 10/060-135
12	20	9	7/2.33	PW 12/020-085	PW 12/020-135
	30			PW 12/030-085	PW 12/030-135
	40			PW 12/040-085	PW 12/040-135
	60			PW 12/060-085	PW 12/060-135