



CHAMPION™ PTCA Guide Wire – Stiff tip for highly stenosed lesions (CTO)



Dimension

Ø 0.014" / 0,36 mm

Length

195 cm
300 cm

Distal tip form

Straight
J-tip

Stent zone stiffness

Standard

Tip stiffness (Gram Force)

3,0gf
4,5gf
6,0gf
9,0gf
12,0gf

Radiopaque length

3 cm

Coil surface

Uncoated – Stainless Steel
Hydrophilic coated

Shaft surface

PTFE coated

Description of item code

PTCA-CTO-S/J-M/HY-3/4,5/6/9/12-D-014-L
(example: PTCA-CTO-S-HY-3-D-014-175)

CTO	For highly stenosed lesions
S	Straight tip
J	J-tip
M	Coil surface uncoated – stainless steel
HY	Coil surface hydrophilic
3	Tip stiffness 3 Gram Force
4,5	Tip stiffness 4,5 Gram Force
6	Tip stiffness 6 Gram Force
9	Tip stiffness 9 Gram Force
12	Tip stiffness 12 Gram Force
D	Stent zone stiffness standard
014:	Diameter of the guide wire 0.014" (0,36 mm)
L:	Length of the guide wire Available in 195, 300 cm