

euca

PW

PTA OTW 0.035" Balloon Catheter

- Dual Lumen Shaft Technology
- Dual Tapered Tip for Powerful Access of Complex Lesions
- Rapid Inflation and Deflation Time



Quality. Made in Germany



euca PW

PTA OTW 0.035" Balloon Catheter

36

Months Shelf Life



Quality. Made in Germany

eucatech AG

Rebgartenweg 27

D-79576 Weil am Rhein / Germany

Phone: +49 7621 4220 - 0

Fax: +49 7621 4220 - 100

E-Mail: eucatech@eucatech.de

Web: www.eucatech.com

eucaPW_EN_0_9_2016 CE0086

Catheter Specifications		
Usable Catheter Length		85 cm / 135 cm
Material	Balloon	Proprietary polyamide compound
	Shaft	Polyamide, dual lumen tube
Shaft Size		5 F - 7 F depending on balloon size
Balloon Folding		Tri-fold
Marker Bands		Two embedded platinum-iridium markers
Tip Entry Profile		0.040 inch (1.02 mm)
Maximum Guide Wire Diameter		0.035 inch (0.89 mm)

Balloon – Ø [mm]	Balloon Length [mm]	Rated Burst Pressure [MPa / atm]	Min. Introducer Sheath [F / mm]	Usable Catheter Length	
				85 cm	135 cm
4	20	1.4 / 14	5 / 1.67	PW04020085	PW04020135
	30			PW04030085	PW04030135
	40			PW04040085	PW04040135
	60			PW04060085	PW04060135
	80			PW04080085	PW04080135
	100			PW04100085	PW04100135
	120			PW04120085	PW04120135
	150			PW04150085	PW04150135
	170			PW04170085	PW04170135
	200			PW04200085	PW04200135
5	20	1.3 / 13	5 / 1.67	PW05020085	PW05020135
	30			PW05030085	PW05030135
	40			PW05040085	PW05040135
	60			PW05060085	PW05060135
	80			PW05080085	PW05080135
	100			PW05100085	PW05100135
	120			PW05120085	PW05120135
	150			PW05150085	PW05150135
	170			PW05170085	PW05170135
	200			PW05200085	PW05200135
6	20	1.2 / 12	5 / 1.67	PW06020085	PW06020135
	30			PW06030085	PW06030135
	40			PW06040085	PW06040135
	60			PW06060085	PW06060135
	80			PW06080085	PW06080135
	100			PW06100085	PW06100135
	120			PW06120085	PW06120135
	150			PW06150085	PW06150135
	170			PW06170085	PW06170135
	200			PW06200085	PW06200135
7	20	1.1 / 11	6 / 2.00	PW07020085	PW07020135
	30			PW07030085	PW07030135
	40			PW07040085	PW07040135
	60			PW07060085	PW07060135
	80			PW07080085	PW07080135
	100		7 / 2.33	PW07100085	PW07100135
	120			PW07120085	PW07120135
	150			PW07150085	PW07150135
	170			PW07170085	PW07170135
	200			PW07200085	PW07200135
8	20	1.1 / 11	6 / 2.00	PW08020085	PW08020135
	30			PW08030085	PW08030135
	40			PW08040085	PW08040135
	60			PW08060085	PW08060135
9	20	1.1 / 11	6 / 2.00	PW08080085	PW08080135
	30			PW09020085	PW09020135
	40			PW09030085	PW09030135
	60			PW09040085	PW09040135
10	20	1.0 / 10	6 / 2.00	PW09060085	PW09060135
	30			PW10020085	PW10020135
	40			PW10030085	PW10030135
	60			PW10040085	PW10040135
12	20	1.0 / 10	7 / 2.33	PW10060085	PW10060135
	30			PW12020085	PW12020135
	40			PW12030085	PW12030135
	60			PW12040085	PW12040135
				PW12060085	PW12060135