

INFLATION DEVICE

OVERVIEW

During minimally invasive balloon dilation procedures, Device can be used to inflate the balloon with fluid, to monitor the pressure in the balloon during the procedure, and to deflate the balloon after dilation.

FEATURES AND BENEFITS

- Quick deflation of every balloon regardless of the size with the 30cc
- Precise pressure increase up to 30ATM
- Registers immediately without pressure decay or drift.

DURABLE LOCKING MECHANISM


- Efficient Intra cam locking mechanism allows to maintain pressure and guarantee fast and automatic release.
- Reliable for multiple inflations

CLEAR POLYCARBONATE SYRINGE

Ensures accurate inflations up to 40ml volume

ANGLED FACE GAUGE

Easy to read for left - or right-handed users



ORDER INFORMATION

CODE NO.	DESCRIPTION	VOLUME (ML)	PRESSURE
F00274	INFLATION DEVICE	40	30 PSI

DISPLAY: Analog
PRESSURE: 30 ATM
STOP COCK: 3-Way 500 PSI (34 Bar)
TUBING: 13" (33.02cm) Braided Extension Tubing with Male Rotating Luer.
BARREL VOLUME: 40ml