


# INTRODUCER SHEATH SET TRACK

**TRACK** is the next generation in vascular access sheath introducers, featuring the latest hemostatic valve technology, for small to large devices.

## STERILE KIT CONTAINS

- One sheath with hemostatic valve and sideport with built-in 3 way stop cock attached to the hub
- One dilator

- One guidewire, L: 50 cm and OD: 0.035" with a flexible straight and J tip and depth markings
- One 18G needle with protector
- One 5-10 cc syringe



## TRACK INTRODUCER SET WITH HEMOSTATIC VALVE FEATURES

- Ergonomically designed hub
- Easy French size identification with color coded hub and 3-way stopcock
- Kink resistant and flexible sheath
- Coated sheath for smooth insertion
- Atraumatic tip of dilator provides smooth transition between dilator and sheath during insertion
- Track version available with a radiopaque atraumatic tip
- Snap fit dilator with french size indicator
- Dilator with curve memory effect for easy sheath placement
- Excellent visibility under fluoroscopy
- Complete hemostasis during insertion and withdrawal of devices

## ORDER INFORMATION

### FEMORAL TYPE

### RADIAL TYPE

CODE NO.	SIZE (F)	LENGTH (CM)	NEEDLE (G)	GUIDE WIRE
F00019	5F	7	20	0.021" / ST - END
F00023	6F	7	20	0.021" / ST - END
F00026	7F	7	20	0.021" / ST - END

CODE NO.	SIZE (F)	LENGTH (CM)	NEEDLE (G)	GUIDE WIRE
F00015	4F	9	18	0.021" / J - END
F00020	5F	11	18	0.035" / J - END
F00024	6F	11	18	0.038" / J - END
F00027	7F	11	18	0.038" / J - END
F00029	8F	11	18	0.038" / J - END
F00031	9F	11	18	0.038" / J - END