


# MANIFOLD

## OVERVIEW

**WELLEX** manifolds help you managing your procedure  
Customized fluid flow control without disconnecting existing lines for connecting new fluid lines.

## ADVANTAGES

Smooth turning handles (ergonomic design)  
Transparent housing for visual control  
Luer lock fitting male/female  
Color coded handles (full-body version)



## FEATURES AND BENEFITS

Manifold product range

**CONTROL VOLUME:** With large bore (.093"/.25cm) inner lumen throughout, including rotator.

**CONTROL ENVIRONMENT** with leak free, smooth turning, low torque handles, with right port orientation and rotating luer.

**CONTROL DEBUBBLING** with clear, polycarbonate bodies.

## CONNECTIONS TO

Catheter  
Saline solution bag  
Contrast medium  
Monitoring systems

## ORDER INFORMATION

CODE NO.	DESCRIPTION	PRESSURE
F00231	2 PORT MANIFOLD RIGTH (ONN/OFF)	500 PSI
F00232	3 PORT MANIFOLD RIGTH (ON/OFF)	500 PSI