

# HEMOSTATIC VALVE SET

## OVERVIEW

The easy push-pull mechanism for both diagnostic & Interventional procedures allows single-handed manipulation to open and close the valve while introducing devices into the vasculature. The hemostasis seal minimizes blood loss while giving the clinician the ability to manipulate interventional devices.

## FEATURES AND BENEFITS

9F compatible lumen allowing delivery of wide range of devices.

Torque and introducer fit 0.014" - 0.021" guide wire.

### PUSH / PULL Y-CONNECTOR

For easy fast and accurate single handed vascular device insertion.  
 Perfect insertion device respecting the integrity of the coated surface.  
 Combination of hemostatic valve and an integrated push-pull valve opener.  
 Clean and safe procedure with a minimized backflow thus limiting blood loss.

### PUSH / CLICK Y-CONNECTOR

Push Click Y-connector provides an easy to use hemostatic device for both diagnostic and interventional procedures.  
 Push Click Y-Connector prevents back-bleeding while balloon catheter/guide wire is inserted into the guiding catheter.  
 Significantly reduces blood loss  
 Push Click Y-connector utilizes a one hand operation both for opening the device by pushing the valve or closing the device by pressing on the lever with a thumb.



## ANGIOPLASTY ACCESSORY COMPONENTS KIT INCLUDES

3 way stopcock  
 Torque device  
 Introducer Needle  
 Y-connector

## ORDER INFORMATION

CODE NO.	DESCRIPTION	COMPATIBLE LUMEN
F00227	PUSH/PULL Y-CONNECTOR WITH SIDE ARM STOPCOCK	9F
F00228	PUSH/PULL Y-CONNECTOR KIT	9F
F00229	PUSH/CLICK Y-CONNECTOR	
F00230	PUSH/CLICK Y-CONNECTOR KIT	